

Air Gel Ankle Stirrup

Material - Polypropylene



Manufactured
in the UK



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

The Air / Gel Stirrup Ankle is comprised of 2 polypropylene outer shells connected with a soft, adjustable foot piece. The air / gel pads can be used for cold therapy to treat oedema and provide uniform compression and a customised fit. The pads are wider at the bottom to prevent biting. It provides stability to the ankle by minimising ankle rotation. Universal size, one size fits all.



Indications for use

- » Acute ankle sprains, chronic ankle instabilities and prophylactic

Features and benefits

- » Provides stability to the ankle by minimising ankle rotation
- » Air/Gel bladder may be use for cold therapy to treat oedema
- » Universal, one size fits all

Sizing Information

Size	Code
Universal	US8813G